




Vermeiren Jazz light transit wheelchair




CRASH TESTED



Lightweight and highly portable



Durable PU wheels



Ergonomic hand grips

The Vermeiren Jazz transit wheelchair is an attendant propelled wheelchair, it has been crash tested and is made to last. Features include ergonomic handgrips, removable armrests and adjustable legrests. Easily folds down for efficient storage and travel.

Swing away footrests (both inward and outwards) and foldable footplates for ease of transfers.

Models

Codes

Vermeiren Jazz light transit wheelchair 460mm wide

464-46

Vermeiren Jazz light transit wheelchair 500mm wide

464-50

SPECIFICATIONS

 1120mm L 935mm H Overall dimensions	 130kg Safe working load
 220mm Armrest height from seat	 420mm Backrest height from seat
 Crash tested Compliance	 min. 370mm max. 440mm Footrest adjustment
 935mm Max height	 16.6kg Product weight
 Belgium Made in	

FEATURES

- Ergonomic handgrips
- Removable armrests
- Lightweight and highly portable
- Swing away footrests (both inward and outwards) and foldable footplates for ease of transfers
- Adjustable and removable legrests
- Attendant propelled folding wheelchair
- 200mm solid PU front wheels and 300mm solid PU rear wheels
- ISO 7176-19 compliant